

Antibodies Product Datasheet

Product Name:	Cdk9 Monoclonal Antibody
Catalog Number:	N0147
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cdk9 expressed in E. Coli.
Concentration:	
Gene Name:	CDK9
Protein Name:	Cell division protein kinase 9
Other Name:	CDK9; CDC2L4; TAK; Cyclin-dependent kinase 9; C-2K; Cell division cycle 2-like protein kinase 4; Cell division protein kinase 9; Serine/threonine-protein kinase PITALRE; Tat-associated kinase complex catalytic subunit
Human Swiss Prot No.:	P50750
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cdk9 Monoclonal Antibody detects endogenous levels of Cdk9 protein.
Storage and Stability:	-20°C/1 year