

Antibodies Product Datasheet

Product Name:	Histone H3 (Tri-Methyl-K10) Antibody
Catalog Number:	D12803
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Methyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Tri-Methyl peptide from human protein at AA range: 10
Concentration:	1 mg/ml
Gene Name:	Histone H3
Protein Name:	Histone H3
Other Name:	H3K10ME3; HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; H
Human Swiss Prot No.:	P68431/Q71DI3/P84243/Q6NXT2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Tri-Methyl-Histone H3 (K10)
Storage and Stability:	-20°C/1 year