

Antibodies Product Datasheet

Product Name:	CYP2D6 Antibody
Catalog Number:	C21046
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CYP2D6.
Concentration:	1 mg/ml
Gene Name:	CYP2D6
Protein Name:	Cytochrome P450 2D6
Other Name:	CYP2D6; CYP2DL1; Cytochrome P450 2D6; CYP11D6; Cytochrome P450-DB1; Debrisoquine 4-hydroxylase
Human Swiss Prot No.:	P10635
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP2D6 Polyclonal Antibody detects endogenous levels of CYP2D6 protein.
Storage and Stability:	-20°C/1 year