

Antibodies Product Datasheet

Product Name:	JAB1 Antibody
Catalog Number:	C11442
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human COPS5. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	COPS5
Protein Name:	COP9 signalosome complex subunit 5
Other Name:	COPS5; CSN5; JAB1; COP9 signalosome complex subunit 5; SGN5; Signalosome subunit 5; Jun activation domain-binding protein 1
Human Swiss Prot No.:	Q92905
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JAB1 Polyclonal Antibody detects endogenous levels of JAB1 protein.
Storage and Stability:	-20°C/1 year