

Antibodies Product Datasheet

Product Name:	UGCGL2 Antibody
Catalog Number:	C50662
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human UGCGL2
Concentration:	1 mg/ml
Gene Name:	UGGT2 UGCGL2 UGT2
Protein Name:	UDP-glucose:glycoprotein glucosyltransferase 2 (UGT2) (hUGT2) (EC 2.4.1.-) (UDP--Glc:glycoprotein glucosyltransferase 2) (UDP-glucose ceramide glucosyltransferase-like 1)
Other Name:	
Human Swiss Prot No.:	Q9NYU1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UGCGL2 at Human
Storage and Stability:	Store at -20° C/1 year