
Antibodies Product Datasheet

Product Name:	S1-Tag Monoclonal Antibody(7B10)
Catalog Number:	N1527
Reactivity:	Species independent
Applications:	WB
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	
Immunogen:	Purified recombinant S1-TAG protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	S1-TAG Monoclonal Antibody detects S1-TAG fusion protein.
Storage and Stability:	-20°C/1 year