

Antibodies Product Datasheet

Product Name:	CALCOCO1 Antibody
Catalog Number:	C49112
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CALCOCO1
Concentration:	1 mg/ml
Gene Name:	CALCOCO1 KIAA1536 PP13275 UNQ2436/PRO4996
Protein Name:	Calcium-binding and coiled-coil domain-containing protein 1 (Calphoglin) (Coiled-coil coactivator protein) (Sarcoma antigen NY-SAR-3)
Other Name:	
Human Swiss Prot No.:	Q9P1Z2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CALCOCO1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year