

Antibodies Product Datasheet

Product Name:	CATW (Cleaved-Val128) Antibody
Catalog Number:	L0417
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CATW (Cleaved-Val128)
Concentration:	1 mg/ml
Gene Name:	CTSW
Protein Name:	CATW (Cleaved-Val128)
Other Name:	Cathepsin W (EC 3.4.22.-) (Lymphopain)
Human Swiss Prot No.:	P56202
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat CATW (Cleaved-Val128, protein was cleaved amino acid sequence between 127-128)
Storage and Stability:	-20°C/1 year