

Antibodies Product Datasheet

Product Name:	Anti-Human Glial Fibrillary Acidic Protein (GFAP) (ABT484) Mouse Monoclonal Antibody
Catalog Number:	V0356
Clone ID:	ABT484
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	GFAP
Gene Name:	P14136
Other Name:	wu:fb34h11, ALXDRD, cb345, etID36982.3, FLJ42474, FLJ45472, GFAP, GFAP_HUMAN, gfapl, Glial fibrillary acidic protein, Intermediate filament protein, wu:fk42c12, xx:af506734, zgc:110485
Isotype	IgG1, Kappa
Localization:	Brain/ Colon
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human GFAP protein.
Storage and Stability:	-20°C/1 year