

Antibodies Product Datasheet

Product Name:	COSMC Antibody
Catalog Number:	C48600
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COSMC
Concentration:	1 mg/ml
Gene Name:	C1GALT1C1 COSMC HSPC067 MSTP143 UNQ273/PRO310
Protein Name:	C1GALT1-specific chaperone 1 (C38H2-like protein 1) (C38H2-L1) (Core 1 beta1,3-galactosyltransferase 2) (C1Gal-T2) (C1GalT2) (Core 1 beta3-Gal-T2) (Core 1 beta3-galactosyltransferase-specific molecular chaperone)
Other Name:	
Human Swiss Prot No.:	Q96EU7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COSMC at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year