

Antibodies Product Datasheet

Product Name:	AP3S2 Antibody
Catalog Number:	C49953
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AP3S2
Concentration:	1 mg/ml
Gene Name:	AP3S2
Protein Name:	AP-3 complex subunit sigma-2 (AP-3 complex subunit sigma-3B) (Adapter-related protein complex 3 sigma-2 subunit) (Sigma-3B-adaptin) (Sigma3B-adaptin) (Sigma-adaptin 3b)
Other Name:	
Human Swiss Prot No.:	P59780
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AP3S2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year