

Antibodies Product Datasheet

Product Name:	Histone H4 (Acetyl-Lys16) Antibody
Catalog Number:	D12251
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p,
Dilutions:	IF: 1:50-200 WB: 1:1000-2000 IHC 1:50-300
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H4 (Acetyl Lys16)
Concentration:	
Gene Name:	HIST1H4A
Protein Name:	Histone H4
Other Name:	HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; HIST2H4A; H4/
Human Swiss Prot No.:	P62805
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H3 (Acetyl Lys16) protein.
Storage and Stability:	-20°C/1 year