

Antibodies Product Datasheet

Product Name:	EpoR (Phospho-Tyr426) Antibody
Catalog Number:	A9036
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human EpoR (phospho Tyr426)
Concentration:	1 mg/ml
Gene Name:	EPOR
Protein Name:	Erythropoietin receptor
Other Name:	EPOR; Erythropoietin receptor; EPO-R
Human Swiss Prot No.:	P19235
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-EpoR (Y426) Polyclonal Antibody detects endogenous levels of EpoR protein only when phosphorylated at Y426.
Storage and Stability:	-20°C/1 year