
Antibodies Product Datasheet

Product Name:	c-Fos Antibody
Catalog Number:	R12-2090
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	c-Fos Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human c-Fos.
Concentration:	1 mg/ml
Gene Name:	FOS
Other Name:	Proto-oncogene c-Fos; Cellular oncogene fos; G0/G1 switch regulatory protein 7; FOS; G0S7
Human Swiss Prot No.:	P01100
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	c-Fos antibody detects endogenous levels of c-Fos.
Storage and Stability:	Stable at -20°C for at least 1 year.