

Antibodies Product Datasheet

Product Name:	TFIIB mouse mAb
Catalog Number:	N0774
Reactivity:	Human, Mouse, Monkey
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human TFIIB protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	gtf2b
Protein Name:	
Other Name:	General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIB; Transcription initiation factor IIB.
Human Swiss Prot No.:	Q00403
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TFIIB and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year