

Antibodies Product Datasheet

Product Name:	MEF-2B Antibody
Catalog Number:	C10843
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEF2B. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	MEF2B
Protein Name:	Myocyte-specific enhancer factor 2B
Other Name:	MEF2B; XMEF2; Myocyte-specific enhancer factor 2B; RSRFR2; Serum response factor-like protein 2
Human Swiss Prot No.:	Q02080
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEF-2B Polyclonal Antibody detects endogenous levels of MEF-2B protein.
Storage and Stability:	-20°C/1 year