
Antibodies Product Datasheet

Product Name:	HDAC7 Antibody
Catalog Number:	R12-2870
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	HDAC7 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HDAC7.
Concentration:	1 mg/ml
Gene Name:	HDAC7
Other Name:	HD7a; HDA7; HDAC7A; Histone deacetylase 7a
Human Swiss Prot No.:	Q8WUI4
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	HDAC7 Antibody antibody detects endogenous levels of HDAC7.
Storage and Stability:	Stable at -20°C for at least 1 year.