

Antibodies Product Datasheet

Product Name:	PTPRO Antibody
Catalog Number:	C51009
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PTPRO
Concentration:	1 mg/ml
Gene Name:	PTPRO GLEPP1 PTPU2
Protein Name:	Receptor-type tyrosine-protein phosphatase O (R-PTP-O) (EC 3.1.3.48) (Glomerular epithelial protein 1) (Protein tyrosine phosphatase U2) (PTP-U2) (PTPase U2)
Other Name:	
Human Swiss Prot No.:	Q16827
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PTPRO at Human, Mouse
Storage and Stability:	Store at -20° C/1 year