

Antibodies Product Datasheet

Product Name:	ARA54 Antibody
Catalog Number:	C15573
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RNF14. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	RNF14
Protein Name:	E3 ubiquitin-protein ligase RNF14
Other Name:	RNF14; ARA54; HRIHFB2038; E3 ubiquitin-protein ligase RNF14; Androgen receptor-associated protein 54; HFB30; RING finger protein 14; Triad2 protein
Human Swiss Prot No.:	Q9UBS8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARA54 Polyclonal Antibody detects endogenous levels of ARA54 protein.
Storage and Stability:	-20°C/1 year