

Antibodies Product Datasheet

Product Name:	FoxO1/3/4 (Phospho-Thr24/32) Antibody
Catalog Number:	A0483
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXO1/3/4-pan around the phosphorylation site of Thr24/32. AA range:15-64
Concentration:	1 mg/ml
Gene Name:	FOXO1/3/4
Protein Name:	Forkhead box protein O1/3/4
Other Name:	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma; FOXO3; FKHRL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1; FOXO4; AFX; AFX1; MLLT7; Forkhead box Q12778/O43524/P98177
Human Swiss Prot No.:	Q12778/O43524/P98177
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-FoxO1/3/4 (T24/32) Polyclonal Antibody detects endogenous levels of FoxO1/3/4 protein only when phosphorylated at T24/32.
Storage and Stability:	-20°C/1 year