

Antibodies Product Datasheet

Product Name:	SKIV2L2 Antibody
Catalog Number:	C48867
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SKIV2L2
Concentration:	1 mg/ml
Gene Name:	SKIV2L2 KIAA0052 Mtr4
Protein Name:	Superkiller viralicidic activity 2-like 2 (EC 3.6.4.13) (ATP-dependent helicase SKIV2L2)
Other Name:	
Human Swiss Prot No.:	P42285
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SKIV2L2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year