

---

## Antibodies Product Datasheet

<b>Product Name:</b>	HRP*Polyclonal Rabbit Anti Human IgM(u chain)
<b>Catalog Number:</b>	SA0101
<b>Reactivity:</b>	HumanIgM(uchain)
<b>Conjugate:</b>	HRP
<b>Applications:</b>	WB,IHC-p,IF/ICC,ELISA
<b>Dilutions:</b>	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
<b>Concentration:</b>	4-5mg/mL
<b>Other Name:</b>	HRP*Polyclonal Rabbit Anti-Human IgM(u chain)
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage and Stability:</b>	-20°C/1 year