

Antibodies Product Datasheet

Product Name:	Phospho-Gab 1 (Y659) Antibody
Catalog Number:	P12-1024
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	Phospho-Gab 1 (Y659) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Gab 1 around the phosphorylation site of Tyr659.
Concentration:	15 mg/ml
Gene Name:	GAB1
Other Name:	GRB2-associated binder 1; GRB2-associated-binding protein 1; growth factor receptor bound protein 2-associated protein 1
Human Swiss Prot No.:	Q13480
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Gab 1 (Phospho-Tyr659) antibody detects endogenous levels of Gab 1 only when phosphorylated at Tyr659.
Storage and Stability:	Store at -20°C/1 year