

Antibodies Product Datasheet

Product Name:	Histone H4 (Acetyl-Lys92) Antibody
Catalog Number:	D12307
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Histone H4 (Acetyl Lys92)
Concentration:	1 mg/ml
Gene Name:	HIST1H4A H4/A H4FA; HIST1H4B H4/I H4FI; HIST1H4C H4/G H4FG; HIST1H4D H4/B H4FB; HIST1H4E H4/J H4FJ; HIST1H4F H4/C H4FC; HIST1H4H H4/H H4FH; HIST1H4I H4/M H4FM; HIST1H4J H4/E H4FE; HIST1H4K H4/D H4FD; HIST1H4L H4/K H4FK; HIST2H4A H4/N H4F2 H4FN HIST2H4; HI
Protein Name:	Histone H4 (Acetyl Lys92)
Other Name:	Histone H4
Human Swiss Prot No.:	P62805
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Histone H4 (Acetyl Lys92)
Storage and Stability:	-20°C/1 year