

Antibodies Product Datasheet

Product Name:	Ron α Antibody
Catalog Number:	C30430
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human MST1R. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	MST1R
Protein Name:	Macrophage-stimulating protein receptor
Other Name:	MST1R; PTK8; RON; Macrophage-stimulating protein receptor; MSP receptor; CDw136; Protein-tyrosine kinase 8; p185-Ron; CD136
Human Swiss Prot No.:	Q04912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ron α Polyclonal Antibody detects endogenous levels of Ron α protein.
Storage and Stability:	-20°C/1 year