

## Antibodies Product Datasheet

<b>Product Name:</b>	Rab2 mouse mAb
<b>Catalog Number:</b>	N0845
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB, ICC
<b>Dilutions:</b>	wb 1:1000 icc 1:200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human Rab2 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	rab2a
<b>Protein Name:</b>	
<b>Other Name:</b>	9330148M11Rik; C80220; GTP binding protein; LHX; OTTMUSP00000004796; RAB2; RAB2 member RAS oncogene family; RAB2A; RAB2A, member RAS oncogene family; RAB2A_HUMAN; Ras related protein Rab 2A; RAS-associated protein RAB2; Ras-related protein Rab-2A; RP23-183L6.1; Small G
<b>Human Swiss Prot No.:</b>	P61019
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Rab2 and does not cross-react with related proteins.
<b>Storage and Stability:</b>	-20°C/1 year