

Antibodies Product Datasheet

Product Name:	Nuclear Receptor NR4A1 (Phospho-Ser351) Antibody
Catalog Number:	A1168sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA(peptide)1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized pospho peptide derived from human Nuclear Receptor NR4A1 (Phospho-Ser351)
Concentration:	1 mg/ml
Gene Name:	NR4A1 GFRP1 HMR NAK1
Protein Name:	Nuclear Receptor NR4A1 (Phospho-Ser351)
Other Name:	Nuclear receptor subfamily 4 group A member 1 (Early response protein NAK1) (Nuclear hormone receptor NUR/77) (Nur77) (Orphan nuclear receptor HMR) (Orphan nuclear receptor TR3) (ST-59) (Testicular receptor 3)
Human Swiss Prot No.:	P22736
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous pospho levels of human Nuclear Receptor NR4A1 (Phospho-Ser351)
Storage and Stability:	-20°C/1 year