
Antibodies Product Datasheet

Product Name:	Desmin Antibody
Catalog Number:	C42160
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Desmin
Concentration:	1 mg/ml
Gene Name:	DES
Protein Name:	Desmin
Other Name:	Desmin; YM33115
Human Swiss Prot No.:	P17661
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.75% sodium azide.
Specificity:	This antibody detects endogenous levels of Desmin at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year