
Antibodies Product Datasheet

Product Name:	Anti-Claudin 3 [ABT-CLD3]
Catalog Number:	V0060
Clone ID:	ABT-CLD3
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	O15551
Gene Name:	CLDN3
Other Name:	Claudin-3; Clostridium perfringens enterotoxin receptor 2; CPE-receptor 2; CPE-R 2; Ventral prostate.1 protein homolog; HRVP1; CLDN3; C7orf1; CPETR2
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Claudin 3 [ABT-CLD3] detects endogenous levels of total Claudin 3 protein.
Storage and Stability:	-20°C/1 year