

Antibodies Product Datasheet

Product Name:	SQSTM1/p62 (Phospho-Ser349) Antibody
Catalog Number:	A9625
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human SQSTM1(Ser349)
Concentration:	1 mg/ml
Gene Name:	SQSTM1 ORCA OSIL
Protein Name:	SQSTM1/p62 (Ser349)
Other Name:	Sequestosome-1 (EBI3-associated protein of 60 kDa) (EBIAP) (p60) (Phosphotyrosine-independent ligand for the Lck SH2 domain of 62 kDa) (Ubiquitin-binding protein p62)
Human Swiss Prot No.:	Q13501
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse SQSTM1/p62 (phospho-Ser349)
Storage and Stability:	-20°C/1 year