

Antibodies Product Datasheet

Product Name:	CtBP1 (Phospho-Ser422) Antibody
Catalog Number:	A8325
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CtBP1 around the phosphorylation site of Ser422. AA range:388-437
Concentration:	1 mg/ml
Gene Name:	CTBP1
Protein Name:	C-terminal-binding protein 1
Other Name:	CTBP1; CTBP; C-terminal-binding protein 1; CtBP1
Human Swiss Prot No.:	Q13363
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-CtBP1 (S422) Polyclonal Antibody detects endogenous levels of CtBP1 protein only when phosphorylated at S422.
Storage and Stability:	-20°C/1 year