

Antibodies Product Datasheet

Product Name:	AKAP 14 Antibody
Catalog Number:	C11279
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKAP14. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	AKAP14
Protein Name:	A-kinase anchor protein 14
Other Name:	AKAP14; AKAP28; A-kinase anchor protein 14; AKAP-14; A-kinase anchor protein 28 kDa; AKAP 28; Protein kinase A-anchoring protein 14; PRKA14
Human Swiss Prot No.:	Q86UN6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AKAP 14 Polyclonal Antibody detects endogenous levels of AKAP 14 protein.
Storage and Stability:	-20°C/1 year