
Antibodies Product Datasheet

Product Name:	CS057 Antibody
Catalog Number:	C42143
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CS057
Concentration:	1 mg/ml
Gene Name:	C19orf57
Protein Name:	Uncharacterized protein C19orf57 (Pre-T/NK cell-associated protein 3B3)
Other Name:	
Human Swiss Prot No.:	Q0VDD7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CS057
Storage and Stability:	-20°C/1 year