

Antibodies Product Datasheet

Product Name:	CIB1 Monoclonal Antibody
Catalog Number:	N0154
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CIB1 expressed in E. Coli.
Concentration:	
Gene Name:	CIB1
Protein Name:	Calcium and integrin-binding protein 1
Other Name:	CIB1; CIB; KIP; PRKDCIP; Calcium and integrin-binding protein 1; CIB; Calcium- and integrin-binding protein; CIBP; Calmyrin; DNA-PKcs- interacting protein; Kinase-interacting protein; KIP; SNK-interacting protein 2-28; SIP2-28
Human Swiss Prot No.:	Q99828
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CIB1 Monoclonal Antibody detects endogenous levels of CIB1 protein.
Storage and Stability:	-20°C/1 year