

Antibodies Product Datasheet

Product Name:	PDZK3 Antibody
Catalog Number:	C10065
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PDZD2. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	PDZD2
Protein Name:	PDZ domain-containing protein 2
Other Name:	PDZD2; AIPC; KIAA0300; PDZK3; PDZ domain-containing protein 2; Activated in prostate cancer protein; PDZ domain-containing protein 3
Human Swiss Prot No.:	O15018
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PDZK3 Polyclonal Antibody detects endogenous levels of PDZK3 protein.
Storage and Stability:	-20°C/1 year