
Antibodies Product Datasheet

Product Name:	VASP (Phospho-Ser157) Antibody
Catalog Number:	A7249
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VASP around the phosphorylation site of Ser157. AA range:124-173
Concentration:	1 mg/ml
Gene Name:	VASP
Protein Name:	Vasodilator-stimulated phosphoprotein
Other Name:	VASP; Vasodilator-stimulated phosphoprotein; VASP
Human Swiss Prot No.:	P50552
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-VASP (S157) Polyclonal Antibody detects endogenous levels of VASP protein only when phosphorylated at S157.
Storage and Stability:	-20°C/1 year