
Antibodies Product Datasheet

Product Name:	Keratin 18 Antibody
Catalog Number:	R12-2950
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Keratin 18 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Keratin 18.
Concentration:	1 mg/ml
Gene Name:	KRT18
Other Name:	CK 18; CK18; CYK18; Cytokeratin endo B; K18; K1C18; KERD; KRT1-18; KRT18; Keratin; Keratin D; Keratin; type I cytoskeletal 18; cytokeratin 18; type I cytoskeletal 18
Human Swiss Prot No.:	P05783
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Keratin 18 Antibody antibody detects endogenous levels of Keratin 18.
Storage and Stability:	Stable at -20°C for at least 1 year.