

Antibodies Product Datasheet

Product Name:	TGF β RIII (Phospho-Thr842) Antibody
Catalog Number:	A1128
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TGF beta Receptor III around the phosphorylation site of Thr842. AA range:801-850
Concentration:	1 mg/ml
Gene Name:	TGFBR3
Protein Name:	Transforming growth factor beta receptor type 3
Other Name:	TGFBR3; Transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGFR-3; Betaglycan; Transforming growth factor beta receptor III; TGF-beta receptor type III
Human Swiss Prot No.:	Q03167
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-TGF β RIII (T842) Polyclonal Antibody detects endogenous levels of TGF β RIII protein only when phosphorylated at T842.
Storage and Stability:	-20°C/1 year