

Antibodies Product Datasheet

Product Name:	BNP Monoclonal Antibody
Catalog Number:	N0078
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Synthetic peptide corresponding to aa (Glu-Pro-Leu-Gln-Glu-Ser-Pro-Arg-Pro-Thr-Gly-Val-Trp-Cys) of human BNP, conjugated to KLH.
Concentration:	
Gene Name:	NPPB
Protein Name:	Natriuretic peptides B
Other Name:	NPPB; Natriuretic peptides B; Gamma-brain natriuretic peptide
Human Swiss Prot No.:	P16860
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	BNP Monoclonal Antibody detects endogenous levels of BNP protein.
Storage and Stability:	-20°C/1 year