

Antibodies Product Datasheet

Product Name:	GTF3C6 Antibody
Catalog Number:	C49960
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GTF3C6
Concentration:	1 mg/ml
Gene Name:	GTF3C6 C6orf51 CDA020 NPD020
Protein Name:	General transcription factor 3C polypeptide 6 (Transcription factor IIIC 35 kDa subunit) (TFIIIC 35 kDa subunit) (TFIIIC35) (Transcription factor IIIC subunit 6)
Other Name:	
Human Swiss Prot No.:	Q969F1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GTF3C6 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year