

Antibodies Product Datasheet

Product Name:	CTXL Antibody
Catalog Number:	C20656
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CTXL. at AA range: 260-290
Concentration:	1 mg/ml
Gene Name:	VSIG2 CTH CTXL UNQ2770/PRO7154
Protein Name:	CTXL
Other Name:	V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein)
Human Swiss Prot No.:	Q96IQ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTXL
Storage and Stability:	-20°C/1 year