

Antibodies Product Datasheet

Product Name:	PPP1A mouse mAb
Catalog Number:	N0833
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PPP1A protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ppp1ca
Protein Name:	
Other Name:	Alpha isoform serine threonine protein phosphatase PP1alpha 1 catalytic subunit; Catalytic subunit; EC 3.1.3.16; MGC15877; MGC1674; PP 1A; PP-1A; PP1A; PP1A_HUMAN; PP1alpha; PP2C ALPHA; PP2CA; Ppp1ca; Protein Phosphatase 2C Alpha Isoform; Serine threonine protein phos
Human Swiss Prot No.:	P62136
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PPP1A and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year