

## Antibodies Product Datasheet

<b>Product Name:</b>	SR-3A Antibody
<b>Catalog Number:</b>	C12016
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human 5-HT-3A. AA range:161-210
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	HTR3A
<b>Protein Name:</b>	5-hydroxytryptamine receptor 3A
<b>Other Name:</b>	HTR3A; 5HT3R; HTR3; 5-hydroxytryptamine receptor 3A; 5-HT3-A; 5-HT3A; 5-hydroxytryptamine receptor 3; 5-HT-3; 5-HT3R; Serotonin receptor 3A; Serotonin-gated ion channel receptor
<b>Human Swiss Prot No.:</b>	P46098
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	SR-3A Polyclonal Antibody detects endogenous levels of SR-3A protein.
<b>Storage and Stability:</b>	-20°C/1 year