

Antibodies Product Datasheet

Product Name:	DOT1L Antibody
Catalog Number:	C43502
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 700-780
Concentration:	1 mg/ml
Gene Name:	DOT1L KIAA1814 KMT4
Protein Name:	Histone-lysine N-methyltransferase, H3 lysine-79 specific (EC 2.1.1.43) (DOT1-like protein) (Histone H3-K79 methyltransferase) (H3-K79-HMTase) (Lysine N-methyltransferase 4)
Other Name:	
Human Swiss Prot No.:	Q8TEK3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DOT1L Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year