

Antibodies Product Datasheet

Product Name:	Claudin 7 (PT0731) mouse mAb
Catalog Number:	N1801
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG2a, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Claudin 7
Concentration:	
Gene Name:	CLDN7 CEPTRL2 CPETRL2
Protein Name:	Claudin-7 (CLDN-7)
Other Name:	
Human Swiss Prot No.:	O95471
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Claudin 7. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year