

Antibodies Product Datasheet

Product Name:	Btk Monoclonal Antibody
Catalog Number:	N0082
Reactivity:	Human,Monkey
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Btk expressed in E. Coli.
Concentration:	
Gene Name:	BTK
Protein Name:	Tyrosine-protein kinase BTK
Other Name:	BTK; AGMX1; ATK; BPK; Tyrosine-protein kinase BTK; Agammaglobulinaemia tyrosine kinase; ATK; B-cell progenitor kinase; BPK; Bruton tyrosine kinase
Human Swiss Prot No.:	Q06187
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Btk Monoclonal Antibody detects endogenous levels of Btk protein.
Storage and Stability:	-20°C/1 year