

Antibodies Product Datasheet

Product Name:	ITK Antibody
Catalog Number:	B1079
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITK
Concentration:	1 mg/ml
Gene Name:	ITK EMT LYK
Protein Name:	ITK
Other Name:	Tyrosine-protein kinase ITK/TSK (EC 2.7.10.2) (Interleukin-2-inducible T-cell kinase) (IL-2-inducible T-cell kinase) (Kinase EMT) (T-cell-specific kinase) (Tyrosine-protein kinase Lyk)
Human Swiss Prot No.:	Q08881
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ITK
Storage and Stability:	-20°C/1 year