

Antibodies Product Datasheet

Product Name:	TPTE Antibody
Catalog Number:	C50762
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TPTE
Concentration:	1 mg/ml
Gene Name:	TPTE
Protein Name:	Putative tyrosine-protein phosphatase TPTE (EC 3.1.3.48) (Cancer/testis antigen 44) (CT44) (Transmembrane phosphatase with tensin homology) (Tumor antigen BJ-HCC-5)
Other Name:	
Human Swiss Prot No.:	P56180
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TPTE at Human
Storage and Stability:	Store at -20° C/1 year