

Antibodies Product Datasheet

Product Name:	EDD Antibody
Catalog Number:	R12-2719
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	EDD Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDD.
Concentration:	1 mg/ml
Gene Name:	UBR5
Other Name:	100 kDa protein; EDD; EDD1; hHYD; Hyperplastic discs protein homolog; KIAA0896; Progesterin-induced protein; Ubiquitin--protein ligase EDD; UBR5
Human Swiss Prot No.:	O95071
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	EDD Antibody antibody detects endogenous levels of EDD.
Storage and Stability:	Stable at -20°C for at least 1 year.