

Antibodies Product Datasheet

Product Name:	VASP (Phospho-Thr278) Antibody
Catalog Number:	A8547
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human VASP (phospho Thr278)
Concentration:	1 mg/ml
Gene Name:	VASP
Protein Name:	Vasodilator-stimulated phosphoprotein
Other Name:	VASP; Vasodilator-stimulated phosphoprotein; VASP
Human Swiss Prot No.:	P50552
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-VASP (T278) Polyclonal Antibody detects endogenous levels of VASP protein only when phosphorylated at T278.
Storage and Stability:	-20°C/1 year